



MORE THAN **1 in 5** AEDs MAY NOT BE READY FOR USE.

Are you *sure* yours are ready?



AEDs are FDA-regulated medical devices approved for use by non-medically trained lay people to treat sudden cardiac arrest. AED owners are responsible for ensuring these life-saving devices are ready to work every time they are needed. But ensuring readiness can be a real challenge for organizations of any size.

Yes, AEDs perform periodic self-tests of internal components, and all AEDs have visual status indicators that alert owners if there's a problem with the device's readiness. But regularly checking for these indications, and keeping track of the various expiration dates of each AED's batteries and pads, are all in human hands. And humans have a lot on their hands already. **But now, there's AED Sentinel, a remote AED monitoring system built for every AED program.**

Today's reliance on human AED inspectors results in a high rate of non-compliance, meaning large numbers of AEDs go uninspected or poorly maintained for long periods of time. This leads to an increased risk of AED failure, avoidable sudden cardiac death and lawsuits. And, of course, there are direct and indirect costs associated with pulling valuable employees away from their primary jobs to perform these checks. *AED Sentinel eliminates the burden and cost of human inspections, and reduces the risk and consequences of AED failure.*

- ➔ **The best automated solution to ensure readiness:** Instantly increase your inspection frequency from (maybe) once per month to multiple checks per day, with 24/7 automated monitoring of your AEDs and accessories.
- ➔ **Cost-effective compliance:** Free up your valuable employees and reduce your AED lawsuit risk with one modest investment. Backed by AED Sentinel Certified Compliance, no matter what brand AEDs you have, you can rest assured you're covered.
- ➔ **True ease of ownership:** Eliminate both the burden and cost of human inspectors, making your AED program truly reliable and giving you peace of mind that your AEDs will be ready every time.
- ➔ **A healthy ROI:** The investment you make in AED Sentinel provides the ultimate return. We deliver on the promise of your AED program by ensuring your devices are always ready when you need them.

About AED Sentinel

AED Sentinel is a technology-based remote AED monitoring system, designed and built by the leading authority in AED program compliance, to provide the ultimate in ease of AED ownership and management. It is thoughtfully designed to provide an uncomplicated solution to a complicated problem, enabling AED program managers to quickly self-install and foolproof their AED programs in minutes. AED Sentinel uses internet-connected optical visual status indicator monitoring sensors, along with an associated software platform and an intuitive user dashboard for monitoring, alerts and reporting. And it works with all current AED makes and models installed in wall-mounted cabinets.

Features and Benefits

- ➔ **Installation:** Easy installation and provisioning allows subscribers to quickly set up and start monitoring.
- ➔ **Alerts:** 24/7 automated monitoring alerts you to:
 - AED cabinet door open
 - AED ready or service required
 - AED batteries and pads within useful life or require replacement
- ➔ **Reporting:** Review, download and/or share program performance reports
 - Access digital documentation of all monitoring results and status check logs in the user-friendly (AED Sentinel Manager)
 - Give Legal/Compliance a copy of your AED Sentinel Certified Compliance certificate, showing how AED Sentinel ensures your AED program is compliant with AED inspection standards, AED laws and manufacturer guidelines
- ➔ **Technology & Integration:**
 - No fixed power source required – hardware uses regular, off-the-shelf batteries that are easily replaced
 - Data transmitted via cellular network, completely independent of your network, so no need to bother the IT team
 - Enables data integration with third-party AED program asset tracking systems
 - AED Sentinel hardware is warranty-protected for the life of your AED Sentinel subscription

Pricing

An AED program you can truly rely on is the ultimate return on investment. AED Sentinel delivers on peace of mind, so you can rest assured that the time and money you've already invested in your AED program has your organization ever-ready. And AED Sentinel eliminates the cost of human inspectors, giving you back the time of your valuable employees. A low annual AED Sentinel subscription fee (with an option to pay monthly) includes hardware and monitoring services so you're all set. And, we back our assurances with *AED Sentinel Certified Compliance* at no additional cost, as well as a hardware warranty for the life of your subscription.

About Readiness Systems

Readiness Systems is your AED program expert. We wrote the national AED Program Design Guidelines that define industry standards, and now we bring technology-based services and solutions to help organizations of all sizes with their AED program readiness and compliance. To learn more, visit www.aedsentinel.com

Contact us or your AED supplier today to ensure your AEDs are *always* ready.

READINESS™
S Y S T E M S

☎ (855) 291-9100
✉ info@aedsentinel.com